



PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____
Comments _____	PO # _____

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE CHAIR

<input type="checkbox"/>	888WAL11-BGRD	JETSTREAM-L	K0823	\$5,495.00
<input type="checkbox"/>	888WAL11-BLUE	JETSTREAM-L	K0823	\$5,495.00

BATTERY

<input type="checkbox"/>	BAT-U1	with seat pan upgrade	K0822	
		12V 35AMP / 50AMP	x2	\$220.00

SEATS (W"xD"xH")

<input type="checkbox"/>	SV201820RB	20x18x20 Recline High-Back Pillow-Top		\$370.00**
<input type="checkbox"/>	SC201822RG	Upgrade 20x18 Recline High-Back Grey Cloth		\$150.00
<input type="checkbox"/>	SP202023RG	Upgrade 20wx20d Seat Pan/Capt back/rcing		\$150.00

Also available in 16" & 18" *Seat Cushion Not Included

ARM REST ASSY

<input type="checkbox"/>	104400-00100-B2	(R) Flipback Adj Arm Rest Assy		\$160.00**
<input type="checkbox"/>	104300-00100-B2	(L) Flipback Adj Arm Rest Assy		\$160.00**

TIRES

<input type="checkbox"/>	TSFR3008	Grey	14" Flat-Free Rear Tires (3.00-8)x2	\$80.00**
--------------------------	----------	------	-------------------------------------	-----------

CASTORS

<input type="checkbox"/>	1351-00601-00	6" Front Castor Assy	x2	\$280.00**
--------------------------	---------------	----------------------	----	------------

ANTI-TIP WHEELS

<input type="checkbox"/>	800405-002-01	Anti-Tip Wheel Assy	x2	Standard
--------------------------	---------------	---------------------	----	----------

FOOT REST

<input type="checkbox"/>	Please Call	Flip Up Foot Plate		Standard
<input type="checkbox"/>	ELRJET-B(K)	Elevating Leg Rest in Blue		\$232.00
<input type="checkbox"/>	ELRJET-R(K)	Elevating Leg Rest in Red		\$232.00
<input type="checkbox"/>	202934-00104 +	(L) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/>	5008-08020-210	Hanger Tightening Knob	x2	\$24.00
<input type="checkbox"/>	202934-00105 +	(R) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/>	5008-08020-210	Hanger Tightening Knob	x2	\$24.00
<input type="checkbox"/>	800907-007-00	Swing Away Foot Rest (R)		\$80.00
<input type="checkbox"/>	800907-007-01	Swing Away Foot Rest (L)		\$80.00

CONTROLLER

<input type="checkbox"/>	1361-00111-00A2	VSI Joystick Controller		\$998.00**
<input type="checkbox"/>	1343-PHFW-SQ	VSI Controller Bracket		\$36.00**

CHARGER

<input type="checkbox"/>	109201-88115	4 amp Charger On-Board		\$230.00**
--------------------------	--------------	------------------------	--	------------

ACCESSORIES

<input type="checkbox"/>	CCH-003(K)	Cane Holder		\$150.00
<input type="checkbox"/>	CH-02(K)	Oxygen Tank Holder (Metal)		\$150.00
<input type="checkbox"/>	WH-0185(K)	Walker Holder		\$150.00
<input type="checkbox"/>	OXY-03(K)	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/>	300909-03	Safety Belt		\$34.00**
<input type="checkbox"/>	109903-00302	Front Castor Curb Guard-Set		\$70.00

TOTAL \$

**STANDARD WITH UNIT

JETSTREAM-L SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....	40"
Base Width.....	26"
Height.....	40"
(Ground to Seat Top)	
Ground Clearance.....	3.5"

WEIGHTS

User Weight Capacity.....	300 Lbs***
Total Unit Weight.....	215 Lbs
Weight of Heaviest Part.....	99 Lbs
Battery Module Weight (U1).....	25 Lbs

TIRES/CASTORS

Front Castor Assy.....	6"
Rear Flat Free.....	14"

BATTERIES

BAT-U1 Module Capacity.....	12v35Ah x2
-----------------------------	------------

CHARGER

Charger-On Board.....	4 amp
-----------------------	-------

MOTOR POWER OUTPUT

Horsepower.....	1.34 x2
-----------------	---------

DRIVING SPECS

Turning Radius.....	26"
Range.....	up to 17.5 miles
Maximum Speed.....	6 mph
Maximum Climbing Angle.....	10% / 6°

***unit may be programmed up to 350 lb wt cap

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.