



888-3L
Sunrunner-3 Deluxe



888-4L
Sunrunner-4 Deluxe

PRICE LIST & ORDER FORM

Effective November 2006

For easy ordering and information contact
Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE SCOOTER

<input type="checkbox"/>	888-3L	SUNRUNNER-3 DELUXE	Burgundy	K0806	\$2,295.00
<input type="checkbox"/>	888-4L	SUNRUNNER-4 DELUXE	Burgundy	K0899	\$2,495.00

BATTERY

<input type="checkbox"/>	BAT-U1	12V 35AMP	x2	\$110.00 ea**
--------------------------	--------	-----------	----	---------------

COMPLIMENTARY BATTERIES ARE NOT WARRANTED

SEATS (W"xD"xH")

<input type="checkbox"/>	SV171/21714RB	171/2"x17" Pillow-Top Recline Seat Vinyl	\$290.00**
<input type="checkbox"/>	SV1717131/2G	Upgrade 17½" Deluxe Captain Seat Vinyl	\$150.00
<input type="checkbox"/>	SC171/21814G	Upgrade 17½" Deluxe Captain Seat Grey Cloth	\$150.00
<input type="checkbox"/>	SV201820RG	Upgrade 20"x18" Recline Capt Grey Vinyl	\$150.00
<input type="checkbox"/>	SP181821RG	Upgrade Pan seat w/captain back recline (18x18)	\$150.00

Note: Pan seat avail in 16x16 & 20x20 / cushion not sold

ARM REST ASSY

<input type="checkbox"/>	104300-88800	(L) Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/>	104400-88800	(R) Arm Rest Assy w/ Arm Pad	\$160.00**

TIRES

<input type="checkbox"/>	TAR26085	Rear	10" Pnuematic Tires-all Models	x2	\$32.00 ea**
<input type="checkbox"/>	205102-88100	Rear	10" Tire Tubes-all Models	x2	\$20.00 ea**
<input type="checkbox"/>	TAF26085	Front	10" Pnuematic Tires-4L	x2	\$32.00 ea**
<input type="checkbox"/>	205102-88100	Front	10" Tire Tubes- 4L	x2	\$20.00 ea**
<input type="checkbox"/>	TAF280250	Front	9" Pnuematic Tires-3L	x1/x2	\$36.00 ea**
<input type="checkbox"/>	205202-77100	Front	9" Tire Tubes-3L	x1/x2	\$20.00 ea**
<input type="checkbox"/>	STU8WS(K)	Front & Rear	888-3L Solid Tire Upgrade		\$200.00
<input type="checkbox"/>	STU8A(K)	Front & Rear	888-4L Solid Tire Upgrade		\$200.00

CHARGER

<input type="checkbox"/>	109201-88114	4 amp Charger On-Board (888-3L & 888-4L)	\$200.00**
--------------------------	--------------	--	------------

ACCESSORIES

<input type="checkbox"/>	209916-88100	Front Wire Mesh-Large Basket	\$50.00**
<input type="checkbox"/>	109901-88101	Front Wicker-Large Basket	\$70.00
<input type="checkbox"/>	CCH-003(K)	Cane Holder	\$150.00
<input type="checkbox"/>	CH-02(K)	Oxygen Tank Holder (Metal)	\$150.00
<input type="checkbox"/>	WH-0185(K)	Walker Holder	\$150.00
<input type="checkbox"/>	WIG-0168	Wig-wag extensions (pair)	\$20.00
<input type="checkbox"/>	RB-0125(K)	Rear Wire Mesh Basket	\$170.00
<input type="checkbox"/>	300909-03	Safety Belt	\$34.00

****STANDARD WITH UNIT**

- All specifications are subject to change without notice.
Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:
- 1) Include the anti-tip wheel or rear castor.
 - 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
 - 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.

TOTAL \$

MIDSIZE SPECIFICATIONS

SUNRUNNER-3 Deluxe(888-3L)

Base LxWxH".....	48-51"x24"x41"
*Ground to Seat Top	
Ground Clearance.....	3"
User Weight Capacity.....	300 Lbs
Total Unit Weight.....	160 Lbs
Battery Module Weight (U1).....	24 Lbs x2
Number, Size Tire (front/rear).....	2.9"/2,10"
BAT-U1 Module Capacity.....	12v35Ah x2
Charger-On Board.....	4 amp
Turning Radius.....	39"
Range.....	up to 24 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% / 8°

SUNRUNNER-4 Deluxe(888-4L)

Base LxWxH".....	49-52"x24"x41"
*Ground to Seat Top	
Ground Clearance.....	3"
User Weight Capacity.....	300 Lbs
Total Unit Weight.....	176 Lbs
Battery Module Weight (U1).....	24 Lbs x2
Number, Size Tire (front/rear).....	2.10"/2,10"
BAT-U1 Module Capacity.....	12v35Ah x2
Charger-On Board.....	4 amp
Turning Radius.....	53"
Range.....	up to 24 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% / 8°